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PRINCIPLE AND INTENDED USE

The VivaChek[™] Ino Blood Glucose Test Strips are thin strips with a chemical reagent system. They work with the VivaChek[™] Ino Blood Glucose Meter as a system to quantitatively measure the glucose concentration in fresh capillary whole blood. The VivaChek[™] Ino blood glucose test is based on the measurement of an electric current generated by the reaction of glucose with the reagents on the electrode of the test strip. The blood sample is drawn into the tip of the test strip through capillary action where glucose in the sample reacts with an enzyme and a mediator. Electrons are generated, producing a current that is proportional to the glucose concentration in the sample. After the reaction time, the glucose concentration in the sample is displayed. The meter is calibrated to display plasma-like concentration results. The VivaChek[™] Ino system is fully compliant with the new EN ISO 15197:2015 International Standard. The VivaChek[™] Ino Blood Glucose Test Strips and VivaChek[™] Ino Blood Glucose Meter are intended for external use only (*in vitro* diagnostic use) for self-testing and health care professional use, as an aid to monitor the effectiveness of diabetes control. The system should not be used for diagnosis of diabetes.

COMPOSITION

Each test strip contains the following reactive chemicals: Glucose oxidase < 25 IU, Mediator < 300 µg. Each test strip vial contains a drying agent. Each test strip pouch contains a drying agent.

STORAGE AND HANDLING

- Store the test strips in a cool, dry place between 2-30°C (36-86°F). Keep away from heat and direct sunlight. Exposure to temperature and / or humidity outside the specified conditions may result in inaccurate readings.
- Do not freeze or refrigerate.
- Use the test strips at temperatures between 5-45°C (41-113°F).
- Use the test strips between 10-90% humidity.
- Do not store the meter, the test strips or control solution near bleach or cleaners that contain bleach.
- Replace the vial cap and close it tightly immediately after removing a test strip.
- Always keep the test strips in the original vial. Use the test strip immediately after removing it from the vial or foil pouch.
- Do not use your test strips beyond the expiry date (printed on the strip vial label or on the foil pouch) or discard date, whichever comes first, because this may cause incorrect test results.
- Note:** All expiry dates are printed in Year-Month format. 2020-01 indicates January, 2020.
- Write the discard date (6 months after first opening the vial) on the vial label when you first open it. Discard any remaining test strips after the discard date.
- Do not use test strips that are torn, bent, or damaged in any way. Do not reuse test strips.
- Keep the vial of test strips away from children. Do not swallow test strips.
- Never ignore symptoms or make significant changes to your diabetes control programme without speaking to your healthcare professional or doctor.

PERFORMING A BLOOD GLUCOSE TEST

Material provided: VivaChek[™] Ino Test Strips
Material required but not provided: VivaChek[™] Ino meter, User Manual, lancing device, clear cap (optional) and a new sterile lancet
Refer to your User Manual for complete instructions for blood sample collection before use.

- Select the puncture site. Wash your hands and the puncture site in warm, soapy water, dry them thoroughly.
 - Prepare the lancing device.
 - Check the expiry date and discard date (printed on the strip vial label or the foil pouch). Do not use test strips beyond the expiry date or discard date.
 - Insert the test strip into the meter in the direction of the arrows. The meter will turn on.
 - Lance the puncture site to obtain a round drop of blood.
 - Touch the blood drop to the strip tip until the meter beeps. Do not apply blood on to the flat surface of the test strip.
 - Your blood glucose test result will appear after the meter counts down from 5 to 1.
- IMPORTANT:** The VivaChek[™] Ino Blood Glucose Monitoring System allows alternative site testing for forearm and palm testing in addition to fingertip testing. There are important differences between forearm, palm and fingertip samples that you should know. Important information about forearm and palm glucose testing:

- When blood glucose levels are changing rapidly such as after a meal, insulin dose or exercise, blood from the fingertips may show these changes more rapidly than blood from other areas.
- Fingertips should be used if testing is within 2 hours of a meal, insulin dose or exercise and any time you feel your glucose levels are changing rapidly.
- You should test with the fingertips anytime there is a concern for hypoglycaemia or you suffer from hypoglycaemia unawareness.

EXPECTED DIABETES CONTROL GOAL

Blood glucose values will vary depending on food intake, medication dosages, health, stress, or exercise. Ideally, you should aim to control your glucose level as close to a normal (non-diabetic) blood glucose level as safely as you can.

The American Diabetes Association suggests the following targets for most non-pregnant adults with diabetes. More or less stringent glycemic goals may be appropriate for each individual. In real life, consult your doctor or health care professional for the target value that is appropriate for you.

Expected blood glucose levels for most non-pregnant adults with diabetes:¹

Time	Range, mg/dL	Range, mmol/L
Before a meal	70 – 130	3.9 – 7.2
1-2 hours after beginning of the meal	Less than 180	Less than 10

Questionable or Inconsistent Results: Repeat the test according to the procedures described in your VivaChek[™] Ino User Manual. If attempts to correct a problem fail, contact your health care professional or doctor.

CHECKING THE SYSTEM

Use only VivaChek[™] Ino Control Solutions.
For complete details about checking the system, refer to your User Manual.
When to check:

- At least once a week
 - When you open a new box of test strips
 - When you want to check the meter and test strips
 - If your test strips were stored in extreme temperature or humidity
 - If you dropped the meter
 - If your test result does not match with how you feel
- For confirmation of results, Control Solution Low tests should fall within the CTRL L range, Control Solution Normal tests should fall within the CTRL N range, and Control Solution High tests should fall within the CTRL H range. When testing with the Control Solution Normal, make sure you are matching the results to the CTRL N range printed on the strip vial label or on the foil pouch or on the box.

CAUTION: If your quality control test result falls outside the control range shown on the test strip vial label or on the foil pouch or on the box, **DO NOT** use the system to test your blood, as the system may not be working properly. If you cannot correct the problem, contact your distributor for help.

LIMITATIONS

- The VivaChek[™] Ino meter, test strips and control solution have been designed, tested and proven to work together effectively to provide accurate blood glucose measurements. Do not use components from other brands.
- Fresh capillary blood may be collected into test tubes containing sodium heparin, lithium heparin if the blood is used within 10 minutes. Do not use sodium fluoride/oxalate or other anticoagulants or preservatives.
- Use only with whole blood. Do not use with serum or plasma samples.
- Very high (above 70%) and very low (below 20%) hematocrit levels can cause false results. Talk to your doctor or health care professional to find out your hematocrit level.
- Abnormally high levels (above 3 mg/dL) of vitamin C and other reducing substances will produce false high blood glucose measurements.
- The system is tested to accurately read the measurement of glucose in whole blood within the range of 10 to 600 mg/dL (0.6-33.3 mmol/L).
- Fatty substances (triglycerides up to 3,000 mg/dL (166.7 mmol/L) or cholesterol up to 500 mg/dL (27.7 mmol/L) have no major effect on blood glucose test results.
- The VivaChek[™] Ino Blood Glucose Monitoring System has been tested and shown to work properly up to 10,000R (3,048 meters).
- Severely ill persons should not run the glucose test with the VivaChek[™] Ino Blood Glucose Monitoring System.
- Patients using oxygen therapy are not recommended for testing with VivaChek[™] Ino Blood Glucose Monitoring System.
- Blood samples from patients in shock, or with severe dehydration or from patients in a hyperosmolar state (with or without ketosis) have not been tested and are not recommended for testing with VivaChek[™] Ino Blood Glucose Monitoring System.
- Dispose of blood samples and materials carefully. Treat all blood samples as if they are infectious materials. Follow proper precautions and obey all local regulations when disposing of materials.

PERFORMANCE CHARACTERISTICS

The VivaChek[™] Ino meter is calibrated by using YSI (Model 2300 STAT PLUS) Glucose Analyzer reference instrument, which is traceable to NIST reference standard.

Repeatability, Precision

Repeatability-Blood		
Interval	Glucose concentration	Standard Deviation (SD) or Coefficient of Variation (CV)
1	40 mg/dL (2.2 mmol/L)	1.5 mg/dL (0.083 mmol/L)
2	80 mg/dL (4.4 mmol/L)	2.6 mg/dL (0.14 mmol/L)
3	130 mg/dL (7.2 mmol/L)	2.5%
4	200 mg/dL (11.1 mmol/L)	3.0%
5	325 mg/dL (18.1 mmol/L)	2.6%
Intermediate Precision-Control Solution		
Interval	Glucose concentration	Standard Deviation (SD) or Coefficient of Variation (CV)
1	40 mg/dL (2.2 mmol/L)	2.0 mg/dL (0.11 mmol/L)
2	120 mg/dL (6.7 mmol/L)	2.8%
3	350 mg/dL (19.4 mmol/L)	2.6%

System Accuracy

The capillary blood glucose measurements from 120 participants were taken by a trained technician using the VivaChek[™] Ino Blood Glucose Meter (y). Capillary blood samples were obtained from fingertip, palm and forearm sampling sites for the VivaChek[™] Ino Blood Glucose Meter testing. Fingertip samples from the same subjects were also analyzed with YSI Model 2300 STAT PLUS Glucose Analyzer (x). The results were compared.

Linear Regression Results: VivaChek [™] Ino (y) vs. YSI Reference (x)				
Sample Site	Slope	Intercept (mg/dL) / (mmol/L)	R	N
Fingertip	0.9847	1.4981/0.0832	0.9969	240
Palm	1.0101	-0.0441/-0.0025	0.9950	202
Forearm	1.0034	-0.0706/-0.0039	0.9935	202

Fingertip samples were used for YSI reference measurement.

The sample range was 23.3 to 534 mg/dL (1.3 to 29.7 mmol/L) for VivaChek[™] Ino Blood Glucose Meter testing with blood sampled from fingertip sites. The sample range was 50.4 to 376 mg/dL (2.8 to 20.9 mmol/L) for VivaChek[™] Ino Blood Glucose Meter testing with blood sampled from palm and forearm sites.

Fingertip Site: System Accuracy Results for Glucose Concentration ≥ 100 mg/dL (5.5 mmol/L)		
Within ± 5%	Within ± 10%	Within ± 15%
142/178 (79.8%)	177/178 (99.4%)	178/178 (100%)
Fingertip Site: System Accuracy Results for Glucose Concentration < 100 mg/dL (5.5 mmol/L)		
Within ± 5 mg/dL (0.28 mmol/L)	Within ± 10 mg/dL (0.56 mmol/L)	Within ± 15 mg/dL (0.83 mmol/L)
56/62 (90.3%)	62/62 (100%)	62/62 (100%)

Palm Site: System Accuracy Results for Glucose Concentration ≥ 100 mg/dL (5.5 mmol/L)		
Within ± 5%	Within ± 10%	Within ± 15%
122/160 (76.3%)	149/160 (93.1%)	159/160 (99.4%)
Palm Site: System Accuracy Results for Glucose Concentration < 100 mg/dL (5.5 mmol/L)		
Within ± 5 mg/dL (0.28 mmol/L)	Within ± 10 mg/dL (0.56 mmol/L)	Within ± 15 mg/dL (0.83 mmol/L)
30/42 (71.4%)	42/42 (100%)	42/42 (100%)

Forearm Site: System Accuracy Results for Glucose Concentration ≥ 100 mg/dL (5.5 mmol/L)		
Within ± 5%	Within ± 10%	Within ± 15%
115/160 (71.9%)	147/160 (91.9%)	158/160 (98.8%)
Forearm Site: System Accuracy Results for Glucose Concentration < 100 mg/dL (5.5 mmol/L)		
Within ± 5 mg/dL (0.28 mmol/L)	Within ± 10 mg/dL (0.56 mmol/L)	Within ± 15 mg/dL (0.83 mmol/L)
31/42 (73.8%)	39/42 (92.9%)	42/42 (100%)

For complete instructions, please refer to the User Manual included with your meter. For any additional questions or issues with this product, please contact your distributor for help.

REFERENCES

- ADA Clinical Practice Recommendations, 2014.

INDEX OF SYMBOLS

	Consult instructions for use		Use by		Contains sufficient for <n> tests
	For <i>in vitro</i> diagnostic use only		Lot number		Control range
	Temperature limitations		Manufacturer		Catalog number
	Do not reuse		Authorized Representative		Use by 6 months from the opening

VivaChek[™]
VivaChek Laboratories, Inc.
913 N Market Street, Suite 200,
Wilmington, DE 19801, USA
www.vivachek.com

EC REP

Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd,
2595AA, The Hague, Netherlands.
E-mail : peter@lotusnl.com
Tel : +3164416899

CE 0197

VivaChek[™] is a trademark of VivaChek Laboratories, Inc.
© 2020 VivaChek Laboratories, Inc.

Number: 1205115601
Effective date: 2020-05-15